
Compassionate use of scuba masks with positive pressure in patients with acute respiratory failure related to COVID-19 infection.



Overview

- Urgent need
- Spanish Red Cross "Channel" Partners (Banco Santander, Airbus, Decathlon)support production and distribution
- Offered as an alternative in RCRC or public system facilities in need



236.259 cases (26th May)



Why?

- Shortage of approved devices
- Development of alternative solutions
- First response in patients with respiratory failure secondary to COVID-19 infection
- Delayed admissions to the ICU
- Accessible/Manageable (Wards, Emergency Room)





Easybreath (Decathlon) Duke (Cressi)

How?

- High O2 Flow + Positive End-Expiratory Pressure (PEEP)+ Scuba mask *
- Autonomous respiratory system (non electrical devices /double oxygen flowmeter)
- PEEP enough to avoid alveolar collapse with improved blood oxygen







How?





Prototype






How?









- Closed system (limit the dispersion of potentially contagious particles)
- Anti-suffocation valve
- Similar materials to authorized masks
- Effect on skin, eyes tested



Initial clinical results

Sample: 22 patients

Age: 63,63 (19 male, 3 female)

ER: fever, dysphnea, Sat 02 <93% (15 l/Min), x-ray COVID-19

Non rebreathing mask.

Sat O2:88,5%.

Blood gas: pH 7.45; p02 53,00; PC02 36,80; HCO3 24,38

Scuba mask.

Sat O2: 96,19% (after 2 mins)

Blood gas: pH 7,46; p02 81,67; pCO2 38,57; HCO3 24,55,(after 1

hour)



Conclusions

- High pressure oxygen system is presented as an effective and safe therapeutic tool.
- Valid alternative in the context of the current exceptional situation in the absence of other approved therapeutic options.
- Low cost production.
- Good results for intermediate care (Spain/effectiveness reported from Italy).
- Homologation in process (regulatory framework).
- Presented to ethics committees.



Humanidad

Imparcialidad

Neutralidad

Independencia

Voluntariado

Unidad

Universalidad



www.cruzroja.es 902 22 22 92